

**NEW!**

# PowerFlex® AFD


## ABSORBENT FOAM DRESSING

**Keeps You in  
the Game!**

- Quick & easy to apply
- High absorption and retention of fluid
- Stops bleeding quickly
- Strong, flexible PowerFlex base
- Won't slip or loosen
- Not affected by sweat
- Covering for MRSA infections
- Effective on lacerations, turf burns and abrasions



### PowerFlex AFD Construction



- ← PowerFlex Layer - Strong, flexible PowerFlex sticks to itself for quick & easy wrapping
- ← Waterproof Film - Prevents leakage and contamination
- ← Absorbent Foam - Foam pad draws in fluid while keeping wound site clean

### TESTIMONIAL

#### Comment after using it on a defensive lineman's sutured hand:

"...a great quality product that eliminates 3 steps in the wound care of lacerations and abrasions. One easy step and the athlete is ready to go. Stays in place with heavy physical activity because Andover Cohesive Bandages do not move. Once the sports medicine community is exposed and trials this new product - this is a product that could possibly skyrocket."

- Neal Hazen  
Head Athletic Trainer  
Ball State University



### ORDERING INFORMATION

U.S. Patent No. 5,762,623 Patent Pending  
U.S. Patent No. 6,566,576 Dicon Technologies

Catalog#	Size	Rolls/Case
3722XX-016	2" x 2.5 yds.	16
3744XX-008	4" x 2.5 yds.	8

Available in (XX): Black (BK) or White (WH).

Sterile (S) = Sterile Product/White only.

3722(S)	2" x 2.5 yds.	16
3744(S)	4" x 2.5 yds.	8



Andover Healthcare, Inc.

9 Fanaras Drive, Salisbury, MA 01952 Toll Free: 800.432.6686 Phone: 978.465.0044 Fax: 978.462.0003

www.andoverhealthcare.com

Proudly Made in USA

**Quick, Easy  
and  
Cost Effective!**

# PowerFlex® AFD

## Absorbent Foam Dressing

U.S. Patent Nos. 5.762.623 Patent Pending U.S. Patent No. 6.566.576 Dicon Technologies

✓ **High  
Absorption**

✓ **Less Dressing  
Changes**

✓ **Easier to  
Apply**

PowerFlex AFD



Traditional Bandaging




**Fluid Retention PowerFlex AFD: 124-133 Grams**

**Fluid Retention 12 Ply Gauze: 3.8-4.3 Grams**

- It would take 29 gauze pads to absorb the same amount of fluid as one PowerFlex AFD bandage.
- PowerFlex AFD expands when wet, increasing absorption and retention of fluid while keeping surrounding tissue healthy. Gauze does not expand when wet.
- PowerFlex AFD has a waterproof layer that prevents leakage and contamination. Traditional bandaging is not waterproof.

### PowerFlex AFD Construction



- ← PowerFlex Layer - Strong, flexible PowerFlex sticks to itself for quick & easy wrapping
- ← Waterproof Film - Prevents leakage and contamination
- ← Absorbent Foam - Foam pad draws in fluid while keeping wound site clean