

PowerFlex®

THE ULTIMATE BANDAGE

NOW AVAILABLE IN CAMOUFLAGE PRINT!

- Strong - 23 lbs tensile strength
- Straight hand tear
- Sweat and water resistant
- Controlled compression - will not constrict
- Quick, easy application - no scissors needed



PowerFlex® ★ CAMO

ORDERING INFORMATION

Catalog#	Size	Rolls/Case
3840CM	4" x 5 yds.	18

U.S. Patent No. 5.762.623 European Patent No. 0839222



Andover Healthcare, Inc.

9 Fanaras Drive, Salisbury, MA 01952

www.andoverhealthcare.com

Toll Free: 800-432-6686

Phone: 978-465-0044

Fax: 978-462-0003

Proudly Made in USA